

# LABEL-FREE CENTRAL MEMORY T CELLS

Isolation and staining with CD8 $\alpha$  and CD62L Fab Streptamers<sup>®</sup>

## INTRODUCTION

T cells play a central role in cell-mediated immunity. The development of cytotoxic T lymphocyte (CTL) responses is necessary for the control of a variety of bacterial and viral infections as well as certain types of malignancies. CTLs, which are largely CD8<sup>+</sup>, traffic to peripheral sites of infection and specifically eliminate their target cells<sup>1</sup>.

In the recent years the therapeutic capacity of these cells has opened up new avenues especially for the development of (personalized) cellular therapies based either on highly enriched primary T cell preparations or purified from *in vitro* cultures. These therapies cover a broad spectrum of applications in order to restore e.g. antiviral immunity or improving strategies to treat certain types of cancer<sup>2-6</sup>.

In order to maintain protection, long-living memory T cells are generated that persist throughout an individual's lifespan<sup>1</sup> without the need to re-

encounter antigen eventually. With the identification of functionally and phenotypically distinct subsets of memory T cells – so called central (T<sub>CM</sub>) and effector memory (T<sub>EM</sub>) T cells – however a division of labor between the different CD8<sup>+</sup> T cell types has become evident<sup>7,8</sup>.

Especially CD8<sup>+</sup>CD45RA<sup>-</sup> T<sub>CM</sub> that are characterized by high expression of the lymph node homing molecule CD62L, have been shown to retain phenotypic as well as functional properties of memory T cells upon adoptive transfer<sup>9-11</sup> and were able to confer protection. Notably this was also observed after transfer of *in vitro* expanded T cell clones that were initially derived from CD62L<sup>+</sup> central memory T cells<sup>10</sup>. Taken together, due to their long-term *in vivo* persistence as well as their conserved functional properties, CD62L<sup>+</sup>CD8<sup>+</sup>CD45RA<sup>-</sup> central memory T cells might be superior for adoptive T cell therapy<sup>12</sup>.






Figure 1. Fab Streptamer<sup>®</sup> complexes for reversible cell isolation (A) or for reversible staining (B).



**Figure 2. Cell isolation with reversible Fab Streptamer® reagents.**

Low-affinity Fab-Streps are reversibly multimerized on Strep-Tactin® microbeads forming a Fab Streptamer® for cell isolation. Treatment of isolated cells with the competing Strep-Tactin® ligand biotin causes disruption of the Fab Streptamer® complex and results in spontaneous dissociation of all monomeric Fab-Streps from the target cell surface.



**Figure 3. Reversible staining with fluorescent CD8 and CD62L Fab Streptamers®.**

CD8 or CD62L fluorescent Fab Streptamers® were used to stain CD8<sup>+</sup> (A) or CD62L<sup>+</sup> cells (B), respectively. Cells were analyzed either before or after treatment with biotin. Remaining Fab-Streps could not be detected after subsequent washing steps using (uncomplexed) PE-labeled Strep-Tactin®. Secondary Fab Streptamer® staining of reversibly stained cells served as control. CD8 and CD62L expression of live CD3<sup>+</sup> T cells is shown. Numbers in dot plots indicate percent of cells within gates.


## RESULTS and DISCUSSION

One difficulty in selecting and isolating the desired target cells is that remaining selection reagents often have a negative impact on the cell product (e.g. by reagent-mediated receptor blockade) often completely eliminating its pivotal function. The Streptamer® technology now circumvents these severe problems, as target cells can be entirely liberated from all components of the Streptamer® reagents. After cell isolation with magnetic Fab Streptamers® (Fig. 1a) or after cell staining with fluorescent Fab Streptamers® (Fig. 1b), Fab Streptamers® can be efficiently removed from the labeled cells by gentle biotin (vitamin H)

mediated dissociation of the Fab-Strep – cell complex from Strep-Tactin® magnetic beads (Fig. 2) or from Strep-Tactin® PE, respectively. Subsequent liberation of the cells from all single low affinity Fab-Streps is then achieved by conventional wash steps, and no remaining Fab-Streps can be detected on the surface of previously stained cells (Fig. 3). In addition, the combined use of reversible Fab Streptamers® directed against CD8 and CD62L now allows for the first time to positively isolate T<sub>CM</sub>s by serial positive magnetic enrichment (Fig. 4). An additional CD45RA depletion of the obtained cell product can also be

used to further discriminate between naïve and memory (CD45RA<sup>-</sup>) T cells. CD62L<sup>+</sup> T cells were highly purified by a double magnetic selection from fresh PBMCs in a two-step approach (Fig. 4). First, CD8<sup>+</sup> cells were positively selected by magnetic CD8 Fab Streptamers<sup>®</sup> and entirely liberated from the reagents. After the second isolation step with magnetic CD62L Fab Streptamers<sup>®</sup>, the CD62L<sup>+</sup>CD8<sup>+</sup> T cell target

population was obtained as a highly pure population. Serial magnetic enrichments can also be performed for other selection strategies. For instance antigen-specific CD8<sup>+</sup> T cells can be selected by the combined use of CD8 Fab Streptamers<sup>®</sup> followed by an antigen-specific selection with MHC I Streptamers<sup>®</sup>.



**Figure 4. Serial magnetic enrichment of CD8<sup>+</sup>CD62L<sup>+</sup> T cells**

Cells were incubated with CD8 magnetic Fab Streptamers<sup>®</sup> for preselection of CD8<sup>+</sup> cells. The resulting positive fraction was then further processed by biotin treatment and subsequent washing to remove all CD8 selection reagents. In a second step, the target CD8<sup>+</sup>CD62L<sup>+</sup> T cell population was then highly enriched from the pre-selected CD8<sup>+</sup> cell pool with CD62L magnetic Fab Streptamers<sup>®</sup>. Live lymphocytes for the respective positive and negative fractions of both selection steps are shown.

## REFERENCES

1. Williams, M.A. & Bevan, M.J. Effector and Memory CTL Differentiation. *Annu Rev Immunol* (2006).
2. Moss, P. & Rickinson, A. Cellular immunotherapy for viral infection after HSC transplantation. *Nat Rev Immunol* 5, 9-20 (2005).
3. Riddell, S.R., et al. Restoration of viral immunity in immunodeficient humans by the adoptive transfer of T cell clones. *Science* 257, 238-241 (1992).
4. Rooney, C.M., et al. Use of gene-modified virus-specific T lymphocytes to control Epstein-Barr-virus-related lymphoproliferation. *Lancet* 345, 9-13 (1995).
5. Rosenberg, S.A., Restifo, N.P., Yang, J.C., Morgan, R.A. & Dudley, M.E. Adoptive cell transfer: a clinical path to effective cancer immunotherapy. *Nat Rev Cancer* 8, 299-308 (2008).
6. Heslop, H.E., et al. Long-term restoration of immunity against Epstein-Barr virus infection by adoptive transfer of gene-modified virus-specific T lymphocytes. *Nat Med* 2, 551-555 (1996).
7. Sallusto, F., et al. Two subsets of memory T lymphocytes with distinct homing potentials and effector functions. *Nature* 401, 708-712 (1999).
8. Sallusto, F., Geginat, J. & Lanzavecchia, A. Central memory and effector memory T cell subsets: function, generation, and maintenance. *Annu Rev Immunol* 22, 745-763 (2004).
9. Turtle, C.J., Swanson, H.M., Fujii, N., Estey, E.H. & Riddell, S.R. A distinct subset of self-renewing human memory CD8+ T cells survives cytotoxic chemotherapy. *Immunity* 31, 834-844 (2009).
10. Berger, C., et al. Adoptive transfer of effector CD8 T cells derived from central memory cells establishes persistent T cell memory in primates. *J Clin Invest* 118, 294-305 (2008).
11. Zhang, Y., Joe, G., Hexner, E., Zhu, J. & Emerson, S.G. Host-reactive CD8+ memory stem cells in graft-versus host disease. *Nat Med* 11, 1299-1305 (2005).
12. Klebanoff, C.A., Gattinoni, L. & Restifo, N.P. CD8+ T-cell memory in tumor immunology and immunotherapy. *Immunity* 21, 214-224 (2006).