

Data Sheet

Strep-Tactin® Superflow®

50% suspension

Cat. No.: 2-1206-025

Version: 14.0
Revision Date: 07.01.2021

Lot No.: 1206-

Description	Immobilized streptavidin variant called Strep-Tactin® (5 mg/ml resin) especially optimized for purification of Strep-tag®II and Twin-Strep-tag® fusion proteins*.
Support	6% agarose, crosslinked
Form	50% suspension in buffer (100 mM Tris/HCl pH 8.0, 1 mM EDTA, 150 mM NaCl, 0.02% sodium azide)
Biotin binding activity	> 300 nmol/ml resin
Stability	6 months after shipping
Storage	recommended: 2—8 °C
Shipping	room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.

Application	To allow an efficient binding of Twin-Strep-tag® or Strep-tag®II proteins to Strep-Tactin® it is strongly recommend using column purification instead of batch purification. It is crucial that protein binding takes place on the column. Even a pre-incubation of resin and lysate prior to filling the resin into a column will lead to decreased protein yields. Batch purification should be performed with MagStrep “type3” XT beads. Further, prolonged batch incubations increase the risk of proteolytic degradation of the target protein including cleavage of the tag.
Elution	1x Buffer E, Strep-Tactin® Elution Buffer with Desthiobiotin: 100 mM Tris/HCl pH 8.0, 150 mM NaCl, 1 mM EDTA, 2.5 mM desthiobiotin. It is possible to use 5—10 mM desthiobiotin to get the target protein eluted at higher concentration.
Regeneration	1x Buffer R, Strep-Tactin® Regeneration Buffer with HABA: 100 mM Tris/Cl, 150 mM NaCl, 1 mM EDTA, 1 mM HABA (hydroxy-azophenyl-benzoic acid). If HABA cannot be efficiently removed from Strep-Tactin® Superflow® by using Buffer W, we recommend using Buffer W at pH 10.5 or alternatively 100 mM Tris base.

* Voss, S. & Skerra, A. (1997) Mutagenesis of a flexible loop in streptavidin leads to higher affinity for the *Strep*-tag II peptide and improved performance in recombinant protein purification. *Protein Eng.* 10, 975-982.

For research use only

Important licensing information

This product is based on Strep-Tactin technology covered by intellectual property (IP) rights and on completion of the sale IBA grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.



Please download all up-to-date manuals, protocols and other material from <http://www.iba-lifesciences.com>.