

Protocol

Classic cloning with pASK-IBA

General protocol for cloning genes into vectors via the multiple cloning site

1 DESCRIPTION

pASK-IBA vectors are devised for *E. coli* expression of recombinant proteins fused with the Strep-tag®II. They include the chloramphenicol resistance cassette for selection of transformed *E. coli* cells, F1 and ColE1 origin for a high plasmid copy number, a multiple cloning site with recognition sites for several common restriction enzymes (e.g., *EcoRI*, *SmaI*, *BamHI*, *SaI*, *PstI*), and the inducible tetracycline promoter/operator for the regulated expression of proteins. Since the *tet*-promoter works independently of the genetic background of *E. coli*, pASK-IBA vectors can be combined with any *E. coli* strain. Expression of the recombinant protein is induced by addition of anhydrotetracycline, which is a derivative of tetracycline that does not exhibit antibiotic activity. Depending on the specific vector, the protein will be localized in either the cytoplasm or periplasm.

2 GENERAL INFORMATION

2.1 PCR with *Pfu* polymerase

PCR instructions given in this protocol are recommended for the use of *Pfu* polymerase. When using another Polymerase than *Pfu*, please refer to the recommendations of the respective manufacturer. "Hot-start" DNA polymerase is inactive until the initial denaturation step of PCR cycling. This reduces non-specific priming or the formation of primer dimers. Using 3'-phosphorothioate (PTO) protected oligonucleotides is recommended to protect against the 3'-exonuclease activity of proof-reading polymerases.

Essential parameters for PCR optimization are the annealing temperature, the duration of synthesis and the template concentration. Primers should have a theoretical melting temperature between 60 °C and 70 °C. Otherwise the primer melting temperatures can be derived by adding the single base melting temperatures of consecutive bases using 4 °C for each GC pairing and 2 °C for each AT pairing (and 1 °C for each GT pairing). Additionally, to the priming sequence a non-priming 5' extension for target vector specific cloning has to be appended to each primer as further specified in Appendix 4.1 for pASK-IBA3C. The annealing temperature should be chosen at least 5 °C below the melting temperature of each primer. The number of cycles should be kept as low as possible to minimize base substitutions.

2.2 Cloning via the multiple cloning site

Foreign genes can be introduced to pASK-IBA vectors after PCR via the multiple cloning site and standard restriction enzymes. The reading frame of the corresponding vector has to be considered if standard restriction enzymes are used. In some vectors with N-terminal Strep-tag[®]II, the tag is followed by the linker sequence 5'-GGCGCC. This sequence is recognized by three different restriction enzymes generating 5'-overhangs. Cleavage with the suitable enzyme and, if necessary, a subsequent filling reaction makes the production of blunt ends possible in all reading frames.

Using standard restriction sites, additional polylinker derived amino acids are appended at the respective end of the recombinant protein. To avoid the fusion of such polylinker derived amino acids pASK-IBA vectors offer a general cloning strategy via Type IIS restriction enzymes, *Bsa*I or *Eco*31I (NEB, Thermo Scientific[™]). They allow the precise fusion of the structural gene with the vector encoded functional elements (depending on the vector, Strep-tag[®]II, OmpA/BM40-signal sequence, protease cleavage site, start codon, or stop codon). To accomplish this, it is necessary to adapt the structural gene at both ends of the coding region via PCR (see cloning scheme at <https://www.iba-lifesciences.com/stargate-cloning/#nav-typeiis-restriction-enzymes>). To avoid the incorporation of base substitutions, PCR should be performed with a proof-reading DNA polymerase (e.g., *Pfu*) using PTO protected primers.

3 PROTOCOL

3.1 PCR with *Pfu* polymerase

The following protocol is based on standard protocols for *Pfu* polymerase PCR.

3.1.1 Mix the following reagents in a PCR tube and a total volume of 50 μ l:

Concentration	Reagent
200 μ M	dNTP
0.1-0.5 μ M	Forward primer
0.1-0.5 μ M	Reverse primer
5 μ M	10x buffer (supplier)
20-200 pg/ μ l plasmid DNA or 0.1-1 ng/ μ l cDNA library	Template DNA
2.5 U	<i>Pfu</i> polymerase*
Ad 50 μ l	Distilled H ₂ O

*Depending on the recommendations of the manufacturer. *Pfu* can also be added after the initial denaturation step

3.1.2 Use a heated lid if available. Alternatively, overlay the sample with mineral oil.

3.1.3 Start temperature cycling:

Temperature	Time	Number of cycles	PCR step
94 °C	180 s	1	Initial denaturation
94 °C	30 s	15-20 cycles for plasmid DNA or 30-40 cycles for cDNA library	Denaturation
55-65 °C	30-60 s		Annealing
72 °C	30-240 s		DNA synthesis
60-72 °C	300 s	1	Terminal extension
4 °C	∞	1	Storage

3.1.4 Quantify the PCR fragment by comparing the band intensity with a DNA standard. Apply two different amounts of PCR product in separate lanes to find a band of equal intensity in the DNA standard, which has to be applied to the same gel.

Alternatively, quantify via photometric measurement.

3.1.5 Isolate the PCR product from the gel. If multiple bands are visible after gel electrophoresis, only isolate the fragment of the expected size. Use of a DNA purification kit to extract the PCR fragment is recommended, since following elution with H₂O restriction with *Eco311* (*Bsa*I) can be performed immediately.

3.2 Cleavage of PCR fragment



The pASK-IBA vectors can be digested with the isoschizomers *Bsa*I or *Eco311*. Since both enzymes show different cutting efficiencies regarding the DNA source (vector DNA or PCR fragment) and incubation time, we determined the cloning efficiency for different digestion protocols. We recommend using *Bsa*I for 1 hour or *Eco311* for 16 hours for the cleavage of both the PCR fragment and the vector. For a detailed cloning scheme see Appendix 4.1.

3.2.1 Mix the following reagents in a PCR tube and a total volume of 50 μ l:

Amount	Reagent
X μ l	PCR fragment in H ₂ O (spin eluate)
5 μ l	10x <i>Eco311</i> (or <i>Bsa</i> I) restriction buffer
10-20 U/ μ g DNA	Restriction enzyme (<i>Eco311</i> or <i>Bsa</i> I)
ad 50 μ l	distilled H ₂ O

3.2.2 Incubate with *Bsa*I at 50 °C for 1 hour (or *Eco311* at 37 °C for 16 hours). Take measures to avoid evaporation.

3.2.3 Purify the desired fragment by using a DNA purification kit.

3.2.4 Quantify PCR fragment concentration to determine the appropriate vector to insert-ratio.

3.3 Cleavage of Vector

3.3.1 Mix the following reagents in a 500 µl PCR tube and a total volume of 50 µl:

Amount	Reagent
2 µg	Vector DNA
5 µl	10x <i>Eco311</i> (or <i>Bsal</i>) restriction buffer
10-20 U	Restriction enzyme (<i>Eco311</i> or <i>Bsal</i>)
ad 50 µl	distilled H ₂ O

3.3.2 Incubate with *Bsal* at 50 °C for 1 hour (or *Eco311* at 37 °C for 16 hours). Take measures to avoid evaporation.

3.3.3 Incubate with *Pst*I for further 30 min at 37 °C or dephosphorylate linearized vector DNA with alkaline phosphatase according to the manufacturer's recommendations. This step is to reduce background colonies which result from re-ligated vector.

3.3.4 Purify vector fragment using a preparative agarose gel and a suitable DNA purification kit.

3.3.5 Quantify vector fragment concentration to determine the appropriate vector to insert-ratio.

3.4 Ligation of PCR fragment and vector



We strongly recommend preparing a negative control without the addition of PCR fragment to quantify background reactions. For ligation use PCR fragment and vector in a molar ratio of 3:1.

3.4.1 Mix the following reagents in a 500 µl reaction tube and a total volume of 20 µl:

Positive control	Negative control	Reagent
100 ng	100 ng	Vector DNA, digested
50 fmol	-	PCR fragment, digested
2 µl	2 µl	Ligation buffer, 10x
1 U	1 U	T4 DNA ligase
ad 20 µl	ad 20 µl	distilled H ₂ O

3.4.2 Incubate overnight at 16 °C. Heat inactivation is not recommended and not necessary.

3.4.3 Store the sample at 4 °C until transformation.

3.4.4 After initial clone selection (DNA mini preparation/restriction analysis), proceed to DNA sequencing.

3.5 Sequencing

For validating correct vector insertion and sequence of the PCR fragment, the following sequencing primers can be used:

Forward: 5' -GAGTTATTTTACCACTCCCT-3'
 Reverse: 5' -CGCAGTAGCGGTAAACG-3'

The sequencing primers are also suitable for cycle sequencing.

4.2 Multiple Cloning sites of pASK-IBA vectors

pASK-IBA2C

```

1      CCATCGAATGGCCAGATGATTAATTCCTAATTTTTGTTGACACTCTATCATTGATAGAGTTATTTTACCACCTCCCTATCA  80
                                     forward primer

                                     OmpA
                                     M K K T A I A
81     GTGATAGAGAAAAGTGAATGAATAGTTCGACAAAAATCTAGATAACGAGGGCAAAAATGAAAAAGACAGCTATCGCGA  160
                                     XbaI

                                     OmpA
I A V A L A G F A T V A Q A G D H G P E F E L G T R G
161    TTGCAGTGGCACTGGCTGGTTTTCGCTACCGTAGCGCAgcccGGAGACCATGGTCCCGAATTCGAGCTCGGTACCCGGGGA  240
                                     BsaI   BsmFI   SstI   KpnI   BamHI
                                     PshAI   EcoRI
                                     NcoI

                                     link   Strep-tag®II
S L E V D L Q G D H G L S A W S H P Q F E K *
241    TCCCTCGAGGTCGACCTGCAGGGGACCATGGTCTCAgcccTTGGAGCCACCCGAGTTCGAAAAATAATAAGCTTGACC  320
      XhoI SalI PstI BsmFI BsaI Eco47III HindIII
      PshAI
      NcoI

321    TGTGAAGTGA AAAATGGCGCACATTGTGCGACATTTTTTTTGTCTGCCGTTTACCGCTACTGCGTCACGGATCTCCACGC  400
                                     reverse primer
    
```

pASK-IBA3C

```

1      CCATCGAATGGCCAGATGATTAATTCCTAATTTTTGTTGACACTCTATCATTGATAGAGTTATTTTACCACCTCCCTATCA  80
                                     forward primer

                                     M G D R G P E
81     GTGATAGAGAAAAGTGAATGAATAGTTCGACAAAAATCTAGATAACGAGGGCAAAAatggGAGACCGCGGTCCCGAAT  160
                                     XbaI
                                     BsaI   BsmFI
                                     PshAI   EcoRI
                                     SacII

                                     link   Strep-tag®II
F E L G T R G S L E V D L Q G D H G L S A W S H P Q F
161    TCGAGCTCGGTACCCGGGATCCCTCGAGGTCGACCTGCAGGGGACCATGGTCTCAgcccTTGGAGCCACCCGAGTTC  240
      SstI KpnI BamHI SalI PstI BsmFI BsaI Eco47III
      SmaI XhoI PshAI
      NcoI

      E K *
241    GAAAAATAATAAGCTTGACCTGTGAAGTGA AAAATGGCGCACATTGTGCGACATTTTTTTTGTCTGCCGTTTACCGCTAC  320
      HindIII reverse primer

321    TCGCTCACGGATCTCCACGCGCCCTGTAGCGGCGCATTAAAGCGCGCGGGTGTGGTGGTTACGCGCAGCGTGACCCTAC  400
    
```

pASK-IBA4C

```

1      CCATCGAATGGCCAGATGATTAATTCCTAATTTTTGTTGACACTCTATCATTGATAGAGTTATTTTACCCTCCCTATCA      80
                                     forward primer

81      GTGATAGAGAAAAGTGAAATGAATAGTTCGACAAAAATCTAGATAACGAGGGCAAAAAATGAAAAAGACAGCTATCGCGA      160
                                     XbaI
                                     M K K T A I A

      OmpA      link      Strep-tag®II      link      R
161      I A V A L A G F A T V A Q A A S W S H P Q F E K G A E
      TTGCAGTGGCACTGGCTGGTTTTCGCTACCGTAGCGCAGGCCGCTAGCTGGAGCCACCCGAGTTCGAAAAAGgcgcCGAG      240
                                     NheI      BbeI      BsaI
                                     EheI      PshAI
                                     KasI
                                     NarI

      D R G P E F E L G T R G S L E V D L Q G D H G L *
      P R S R I R A R Y P G I P R G R P A G G P W S L I S
      T A V P N S S S V P G D P S R S T C R G T M V S D I *
241      ACCGCGGTCCCGAATTCGAGCTCGGTACCCGGGATCCCTCGAGGTCGACCTGCAGGGGGACCATGGTCTCTgataTCTA      320
      SacII      EcoRI      KpnI      BamHI      SalI      PstI      BsmFI      BsaI      EcoRV
      BsmFI      SstI      SmaI      XhoI      PshAI
                                     NcoI

      N *

321      ACTAAGCTTGACCTGTGAAGTGAAAAATGGCGCACATTGTGCGACATTTTTTTGTCTGCCGTTTACCGCTACTGCGTCA      400
      HindIII      reverse primer
    
```

pASK-IBA5C

```

1      CCATCGAATGGCCAGATGATTAATTCCTAATTTTTGTTGACACTCTATCATTGATAGAGTTATTTTACCCTCCCTATCA      80
                                     forward primer

81      GTGATAGAGAAAAGTGAAATGAATAGTTCGACAAAAATCTAGATAACGAGGGCAAAAAATGGCTAGCTGGAGCCACCCGC      160
                                     XbaI      link      Strep-tag®II
                                     M A S W S H P
                                     NheI

      D R G P E F E L G T R G S L E V D L Q G
      link R P R S R I R A R Y P G I P R G R P A G G
161      Q F E K G A E T A V P N S S S V P G D P S R S T C R G
      AGTTCGAAAAAGgcgcCGAGACCCGGTCCCGAATTCGAGCTCGGTACCCGGGATCCCTCGAGGTCGACCTGCAGGGGG      240
      BbeI      BsaI      BsmFI      SstI      KpnI      BamHI      SalI      PstI      BsmFI
      EheI      PshAI      EcoRI      SmaI      XhoI      PshAI
      KasI      SacII
      NarI

      D H G L *
      P W S L I S N *
      T M V S D I *
241      ACCATGGTCTCTgataTCTAACTAAGCTTGACCTGTGAAGTGAAAAATGGCGCACATTGTGCGACATTTTTTTGTCTGC      320
      NcoI      EcoRV      HindIII
      BsaI

321      CGTTTACCGCTACTGCGTCACGGATCTCCACGCGCCCTGTAGCGGCGCATTAAAGCGCGGGGTGTGGTGGTTACGCGCA      400
      reverse primer
    
```

pASK-IBA6C

```

1      CCATCGAATGGCCAGATGATTAATTCCTAATTTTTGTTGACACTCTATCATTGATAGAGTTATTTTACCACCTCCCTAT 78
                                     forward primer

                                     M K K T A I A
79     CAGTGATAGAGAAAAGTGAAATGAATAGTTCGACAAAAATCTAGATAACGAGGGCAAAAAATGAAAAAGACAGCTATCGC 158
                                     XbaI

                                     OmpA          link          Strep-tag®II          Factor Xa
159     I A V A L A G F A T V A Q A A S W S H P Q F E K I E 238
GATTGCAGTGGCACTGGCTGGTTTCGCCTACCGTAGCGCAGGCCGCTAGCTGGAGCCACCCGAGTTCGAAAAAATCGAAG
                                     NheI

                                     R P R S R I R A R Y P G I P R G R P A G G P W S
                                     E T A V P N S S S V P G D P S R S T C R G T M V S
239     G R R D R G P E F E L G T R G S L E V D L Q G D H G L 318
GgcgcCGAGACCGGGTCCCGAATTCGAGCTCGGTACCCGGGATCCCTCGAGGTCGACCTGCAGGGGACCATGGTCTC
BbeI BsaI BsmFI SstI KpnI BamHI Sali PstI BsmFI BsaI
EheI PshAI EcoRI SmaI XhoI PshAI
KasI SacII NcoI
NarI

L I S N *
  D I *
  *
319     TgataTCTAACTAAGCTTGACCTGTGAAGTGAAAAATGGCGCACATTGTGCGACATTTTTTTTGTCTGCCGTTTACCGCT 398
EcoRV HindIII reverse primer

399     ACTGCGTCACGGATCTCCACGCGCCCTGTAGCGCGCATTAAGCGCGGGTGTGGTGGTTACGCGCAGCGTGACCGCT 478

```

pASK-IBA7C

```

1      CCATCGAATGGCCAGATGATTAATTCCTAATTTTTGTTGACACTCTATCATTGATAGAGTTATTTTACCACCTCCCTATCA 80
                                     forward primer

                                     link Strep-tag®II
79     GTGATAGAGAAAAGTGAAATGAATAGTTCGACAAAAATCTAGATAACGAGGGCAAAAAATGGCTAGCTGGAGCCACCCGC 160
                                     XbaI                                     NheI

                                     R P R S R I R A R Y P G I P R G R P
                                     factor Xa          E T A V P N S S S V P G D P S R S T C
161     Q F E K I E G R R D R G P E F E L G T R G S L E V D L 240
AGTTCGAAAAAATCGAAGGgcgcCGAGACCGGGTCCCGAATTCGAGCTCGGTACCCGGGATCCCTCGAGGTGACCTG
BbeI BsaI BsmFI SstI KpnI BamHI Sali PstI
EheI PshAI EcoRI SmaI XhoI
KasI SacII
NarI

A G G P W S L I S N *
  R G T M V S D I *
  Q G D H G L *
241     CAGGGGACCATGGTCTCTgataTCTAACTAAGCTTGACCTGTGAAGTGAAAAATGGCGCACATTGTGCGACATTTTTTT 320
BsmFI BsaI EcoRV HindIII
PshAI
NcoI

321     TGTCTGCCGTTTACCGCTACTGCGTACCGATCTCCACGCGCCCTGTAGCGCGCATTAAGCGCGGGGTGTGGTGGTT 400
                                     reverse primer

```




Check our Downloads page

www.iba-lifesciences.com/download-area.html

for the latest version of this manual.



Info on warranty / licensing and trademarks available at:

www.iba-lifesciences.com/patents-licenses-trademarks.html



If you have any questions, please contact

strep-tag@iba-lifesciences.com

We are here to help!

