


Strep-Tactin® HRP

Cat. No.: 2-1502-001, 2-1502-002

Description	Strep-Tactin® labelled with horse radish peroxidase (HRP, 2-2.5 moles per mole Strep-Tactin®). Strep-Tactin® is a streptavidin variant specially optimized for Strep-tag®II and Twin-Strep-tag® fusion protein binding.
Form	Suspension in PBS containing approximately 5 mg/ml bovine serum albumin as well as stabilizing agents (see MSDS for details).
Peroxidase activity	800-1100 ABTS units per mg conjugate (Unit definition: one unit is the enzyme activity which oxidizes 1 µmole 2,2'-azino-bis(3-azino-benz-thiazoline-6-sulfonic acid) per min under assay conditions (25 °C, pH 5)
Amount	500 µl (2-1502-001), 20 µl (2-1502-002)
Concentration	0.95-1.05 mg/ml conjugate (Bradford dye binding assay)
Stability	12 months after shipping
Storage	Recommended: 2-8 °C
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Application	<p>Western Blot: dilute 1:4 000 or higher for chromogenic detection and 1:100 000 for chemiluminescence detection</p> <p>ELISA: dilute between 1:2 000 and 1:10 000 for chromogenic detection; start with 1:5 000.</p> <div style="display: flex; align-items: flex-start;">  <div style="margin-left: 20px;"> <p>The detection limit may vary depending on protein size. We recommend starting with a concentration that will yield an absorbance of roughly 0.5. The data provided on the left was obtained with Twin-Strep-tag® fused proteins in an ELISA, with a Strep-Tactin® HRP dilution of 1:5 000.</p> </div> </div>
Manufacturer	IBA Lifesciences

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.