

**! For research use only**

Protocol

# CD19 Nano-TACS<sup>®</sup> Agarose Column Starter Kit

Cat. no. 6-3307-002

mouse, for splenocytes

## 1. GENERAL INFORMATION & TECHNICAL SPECIFICATIONS

Kit components:

Cat. no.	Product	Quantity	Required/isolation
6-6310-300	Strep-Tactin <sup>®</sup> TACS Agarose Column, 0.3 ml	2	1
6-8507-150	CD19 Nano-Strep, mouse, lyophilized, 50 µg	1	20 µg
6-6996-001	Biotin stock solution, 100 mM, 250 µl	1	60 µl
6-6320-025	10x Buffer CI, 25 ml 10x PBS containing 10 mM EDTA and 5% BSA	1	~2 ml
6-3333-001	TACS Column Adapter (0.3 ml column)	1	1

**Required:** ddH<sub>2</sub>O for Buffer CI dilution

**Column specifications:** **2 x 10<sup>7</sup> target cells** out of 1 x 10<sup>8</sup> splenocytes  
Reservoir volume: 3 ml; **For single use only!**

**Storage:** Store all components at 2 - 8 °C. Store reconstituted Nano-Strep at -80 °C.  
(Buffer CI may also be stored at 15 - 25 °C)

**Stability:** 6 months after shipping.

**Shipping:** Room temperature

**Hazards:** Products are not classified as hazardous according to (EC) No 1272/2008 [CLP].  
Material Safety Data Sheets are provided.

## 2. INITIAL PREPARATIONS

Allow the reagents and the sample to equilibrate to room temperature (RT) prior to use. For a sterile isolation, work under a safety cabinet. **Reagent volumes will be sufficient for one selection process.**

### 2.1. Reagent preparation

- 2.1.1.** Prepare 1x Buffer CI from 10x stock by diluting with ddH<sub>2</sub>O. Degas buffer before use, as air bubbles could block the column.
- 2.1.2.** Dissolve **one vial** of lyophilized Fab-Strep (**50 µg**) in **1 ml** Buffer CI by carefully pipetting up and down (avoid foam formation). **Do not vortex!**



Required per column: **20 µg** Fab-Strep in **400 µl** Buffer CI. Store remaining Fab-Strep solution at -80 °C (stable for 6 months) if not required immediately

- 2.1.3.** Prepare 1 mM Biotin Elution Buffer by adding **60 µl** of the 100 mM Biotin stock solution to **6 ml** Buffer CI. Mix thoroughly.

### 2.2. Sample preparation


- 2.2.1.** Prepare splenocytes in Buffer CI.
- 2.2.2.** Adjust the concentration of the splenocytes to 1 – 5 x 10<sup>7</sup> total cells per ml. To remove clumps and to prevent aggregates, pass splenocytes through a 40 µm nylon mesh before isolation.

### 2.3. Column preparation



**2.3.1. Remove** the caps at the top and at the bottom of the column. Allow the storage solution to drain. Place the Strep-Tactin® TACS Agarose Column into the TACS Column Adapter.

**2.3.2. Wash** the Strep-Tactin® TACS Agarose Column by applying **1 ml** Buffer CI and allow the buffer solution to enter the packed bed completely.



**2.3.3. Load** the **400 µl** Nano-Strep solution (2.1.2.) onto the Strep-Tactin® TACS Agarose Column. Let the Nano-Strep solution enter the packed bed completely. Incubate for **2 min**.


**2.3.4. Wash** the Strep-Tactin® TACS Agarose Column with **600 µl** Buffer CI. Discard effluent and change collection tube. Strep-Tactin® TACS Agarose Column is now ready for cell isolation.



**Do not interrupt the procedure for more than 60 min.**

## 3. PROTOCOL

### 3.1. Cell isolation from splenocytes




**3.1.1. Load**  
Apply diluted splenocytes (2.2.2.) in steps of **max. 3 ml**. Collect flow-through containing unlabeled cells.



### 3.1.2. Wash

Apply **3x 3 ml** Buffer CI. (In each step: Let the buffer solution enter the gel bed completely). The agarose bed should now be white again.



### 3.1.3. Elute

From this step on your effluent contains your target cells. Use a **new collection tube**. Apply **400 µl** Biotin Elution Buffer (2.1.3.) and incubate for **5 min**. Elute target cells by applying **2.6 ml** Biotin Elution Buffer. Elute a second time with additional **3 ml** Biotin Elution Buffer.

## 3.3. Further procedure

Centrifuge your eluted cell suspension for **10 min** at **300 x g**. Discard the supernatant and dissolve cell pellet in your desired buffer.



If you plan to continue with a biotin-sensitive assay, please remove biotin by washing with **50 ml** Buffer CI twice. Discard supernatant **completely**.

## 4. TROUBLESHOOTING

### Low yield

#### Option 1:

Check for biotin contamination in your samples.

#### Option 2:

Use flow restrictor during sample loading.

#### Option 3:

Re-apply flow-through (depleted sample) to the column (3.1.1./3.2.1.).


### Low purity

Invert columns after each wash step three times.


## 5. EXAMPLE DATA

Separation of CD19<sup>+</sup> B cells from mouse spleen using the CD19 Nano-TACS® Agarose Column Starter Kit. Unlysed cells were stained with anti-mouse CD19-FITC (1D3) / CD45-eF450 and analyzed by flow cytometry (CytoFlex, BC). Dead cells were excluded from the analysis using DAPI staining. Doublet and debris discrimination were performed using FSC/SSC signals.

### Before isolation



### After isolation





Watch this How-to video to see an exemplary isolation

<https://www.youtube.com/watch?v=bAaMgNikTDI>



Check our Downloads page

<https://www.iba-lifesciences.com/resources/download-area/>

for the latest version of this protocol



Info on warranty / licensing and trademarks available at:

[www.iba-lifesciences.com/patents-licenses-trademarks/](http://www.iba-lifesciences.com/patents-licenses-trademarks/)



If you have any questions, please contact

[strep-tag@iba-lifesciences.com](mailto:strep-tag@iba-lifesciences.com)

We are here to help!

