



## Compatible reagents for Strep-Tactin®

Indicated detergent concentrations refer to the actual tested values (does not correspond to the maximum value).

Reagent	Concentration	Affinity tag
<b>Reducing Agents</b>		
DTT	50 mM	Strep-tag
β-Mercaptoethanol	50 mM	Strep-tag
TCEP	10 mM	Strep-tag
<b>Non-Ionic Detergents</b>		
Octyltetraoxyethylene	C8E4	max. 0.88 %
Decylpentaoxyethylene	C10E5	0.12 %
APO106 - Anapoe	C10E6	0.03 %
Octaethylene glycol monododecyl ether	C12E8	0.005 %
Dodecyl nonaoxyethylene (Thesit)	C12E9	0.023 %
Decyl-β-D-maltoside	DM	0.35 %
N-dodecyl-β-D-maltoside	DDM; LM	0.1 %
Digitonin		0.50 %
N-nonyl-β-D-glucopyranoside	NG	0.2 %
Nonidet P40	NP40	2 %
N-octyl-β-D-glucopyranoside	OG	2.34 %
Triton X-100	TX	2 %
TritonX-114		2 %
Tween20		2 %
Igepal CA-630		0.5 %
Brij-35		0.1 %
<b>Ionic Detergents</b>		
N-lauryl-sarcosine	2 %	Strep-tag
8-HESO; N-octyl-2-hydroxy-ethylsulfoxide	1.32 %	Strep-tag
Sodium-N-dodecyl sulfate	SDS	0.1 %
Sodium deoxycholate (RIPA)		0.5 %
		Twin-Strep

<b>Zwitter-Ionic Detergents</b>			
3-[(3-cholamidopropyl) dimethylammonio]-1-propanesulfonate	CHAPS	0.5 %	Strep-tag
N-decyl-N,N-dimethylamine-N-oxide	DDAO	0.034 %	Strep-tag
Lauryldimethylamine-oxide	LDAO	0.13 %	Strep-tag
Fos-Choline-12		0.1 %	Strep-tag
<b>Others</b>			
Ammonium sulfate	(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	2 M	Strep-tag
CaCl <sub>2</sub>		max. 1 M	Strep-tag
EDTA		50 mM	Strep-tag
EGTA		5 mM	Strep-tag
Ethanol		10 %	Strep-tag
Guanidine		max. 1 M	Strep-tag
Glycerol		max. 25 %	Strep-tag
Imidazole		max. 250 mM	Strep-tag
MgCl <sub>2</sub>		1 M	Strep-tag
NaCl		5 M	Strep-tag
KCl		25 mM	Strep-tag
Urea		max. 1 M	Strep-tag
Cholesterol hemisuccinate	CHS	0.1 %	Strep-tag
<b>Buffers</b>			
HEPES			Strep-tag
Phosphate buffered saline	PBS		Strep-tag
PIPES-NaOH, pH 7.2		10 mM	Twin-Strep
Potassium phosphate buffer		50 mM	Strep-tag
Radioimmunoprecipitation assay buffer	RIPA		Twin-Strep
Tris-HCl		100 mM	Strep-tag