



IBA Lifesciences GmbH
Rudolf-Wissell-Str. 28
37079 Goettingen
Germany
Tel.: +49 (0) 551-5 06 72-0
E-mail: info@iba-lifesciences.com
www.iba-lifesciences.com

Data Sheet

Strep-Tactin® TACS Agarose Column

Cat. No.: 6-6310-300, 6-6310-001

Version: 2.1
Revision Date: 28.02.2024

Description	Ready-to-use column containing Strep-Tactin® (1-2 mg/ml resin)-coated agarose developed for the isolation of cells (Traceless affinity cell selection – TACS).
Support	Cross-linked agarose
Form	Pre-packed in buffer: Phosphate buffered saline (PBS), 1 mM EDTA, 0.01% sodium azide.
Biotin binding activity	> 40 nmol/ml resin
Stability	6 months after shipping
Storage	Recommended: 2-8 °C
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Application	Strep-Tactin® TACS Agarose column is suitable for <ul style="list-style-type: none">- Direct isolation of cells expressing a surface protein fused to a Strep-tag®II or Twin-Strep-tag®- Cell isolation mediated by a Strep-tag®II- or Twin-Strep-tag®-fusion protein- Negative cell selection or depletion (e.g. using biotinylated antibodies)
Recommended usage	Cell isolation via a Twin-Strep-tag®-fusion protein: load 50 µg of a 50 kDa protein per 1 ml bed volume. Elute with 1 mM biotin. Columns are for single use only.

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.