

# Cell isolation with Strep-Tactin® TACS Agarose columns (0.3 ml)

for splenocytes or other single cell suspensions

## 1. Required material

| Cat. No.   | Product  | Required/isolation |
|------------|--|--------------------|
| 6-6310-300 | Strep-Tactin® TACS Agarose Column, 0.3 ml                                | 1                  |
| -          | Protein of choice fused to a Strep-tag®II or Twin-Strep-tag®             | 20 µg*             |
| 6-6325-001 | Biotin stock solution, 100 mM, 1 ml                                      | 60 µl              |
| -          | Buffer for cell isolation, e.g. 1x PBS containing 1 mM EDTA and 0.5% BSA | ~ 20 ml            |
| 6-3333-001 | TACS Column Adapter (1 ml column)  | 1                  |

\*Amounts refer to a 50 kDa Twin-Strep-tag® protein

## 2. Initial preparations

### 2.1. Reagent preparation

Allow the reagents to equilibrate to **room temperature** (RT) prior to use. For a sterile isolation, work under a safety cabinet. **The following volumes will be sufficient for one selection process.**



Watch this How-to video to see an exemplary isolation:

<https://www.youtube.com/watch?v=bAaMgNikTDI>



- 2.1.1.** Prepare a cell isolation buffer of choice, e.g. 1x PBS containing 1 mM EDTA and 0.5% BSA. The buffer should not contain biotin. Degas buffer before use, as air bubbles could block the column.
- 2.1.2.** Dilute your Strep-tag®II or Twin-Strep-tag® fusion protein to a concentration of **50 µg/ml** with cell isolation buffer.



**Titration** of optimal isolation conditions might be necessary. The following instructions are an example for isolating cells with a **50 kDa** protein fused to a Twin-Strep-tag®.

- 2.1.3.** Prepare **1 mM** Biotin Elution Buffer by diluting **60 µl** of 100 mM Biotin stock solution in **6 ml** cell isolation buffer. Mix thoroughly.

### 2.2. Sample preparation

- 2.2.1.** Prepare splenocytes in cell isolation buffer.
- 2.2.2.** Adjust the concentration of the splenocytes to up to **5 x 10<sup>7</sup> total cells per ml**. To remove clumps and to prevent aggregates, pass splenocytes through a 40 µm nylon mesh before isolation.

### 2.3. Column preparation



**2.3.1. Remove** the cap and cut the sealed end of the column at notch. Allow the storage solution to drain. Place the Strep-Tactin® TACS Agarose Column into the TACS Column Adapter.



**2.3.2. Wash** the Strep-Tactin® TACS Agarose Column by applying **1 ml** cell isolation buffer and allow the buffer solution to enter the packed bed completely.



**2.3.3. Load 400 µl** Twin-Strep-tag® fusion protein (2.1.2.) onto the Strep-Tactin® TACS Agarose Column. Let the protein solution enter the packed bed completely. Incubate for **2 min**.




**2.3.4. Wash** the Strep-Tactin® TACS Agarose Column with **600 µl** cell isolation buffer. Discard flow-through and change collection tube. The Strep-Tactin® TACS Agarose Column is now ready for cell isolation.



Do not interrupt the procedure for **more than 60 min**.

## 3. Protocol

### 3.1. Cell isolation from splenocytes



#### 3.1.1. Load

Apply diluted splenocytes (2.2.2.) in steps of **max. 3 ml**. Collect flow-through containing unlabeled cells.



#### 3.1.2. Wash

Apply **3x 3 ml** cell isolation buffer. (In each step: Let the buffer solution enter the gel bed completely). The agarose bed should now be white again.



#### 3.1.3. Elute

From this step on your flow-through contains your target cells. Use a new collection tube. Apply **400 µl** Biotin Elution Buffer (2.1.3.) and incubate for **5 min**. Elute target cells by applying **2.6 ml** Biotin Elution Buffer. Elute a second time with additional **3 ml** Biotin Elution Buffer.

### 3.2. Further procedure

Centrifuge your eluted cell suspension for **10 min** at **300 x g**. Discard the supernatant and dissolve cell pellet in your desired buffer.



If you plan to continue with a biotin-sensitive assay, please remove biotin by washing with **50 ml** cell isolation buffer twice. Discard supernatant completely.

## 4. Troubleshooting

### Low yield

|                 |  |
|-----------------|--|
| <i>Option 1</i> | Check for biotin contamination in your samples.        |
| <i>Option 2</i> | Re-apply flow-through (depleted sample) to the column. |


### Low purity

Invert columns after each wash step three times.

## 5. Example data

Separation of CD3<sup>+</sup> T cells from mouse spleen using a Fab Fragment against CD3 fused to a Twin-Strep-tag<sup>®</sup>. Unlysed cells were stained with anti-mouse CD3-PE (145-2C11) / CD45-eF450 (1D4) and analyzed by flow cytometry (CytoFlex, BC). Dead cells were excluded from the analysis using DAPI staining. Doublet and debris discrimination were performed using FSC/SSC signals.

### Before isolation



### After isolation

