

Data Sheet

Competent *E. coli* Top10

Cat. No.: 5-1600-020

Version: 8.0
 Revision Date: 03.06.2021

Description	Competent <i>E. coli</i> TOP10 cells are ready for heat shock transformation with vector DNA and its subsequent propagation for cloning and transfection purposes.
Genotype	F ⁻ <i>mcrA</i> Δ(<i>mrr-hsdRMS-mcrBC</i>) φ80 <i>lacZ</i> Δ <i>M15</i> Δ <i>lacX74</i> <i>recA1</i> <i>ara</i> Δ <i>139</i> Δ(<i>ara-leu</i>)7697 <i>galU</i> <i>galK</i> <i>rpsL</i> (Str ^R) <i>endA1</i> <i>nupG</i>
Transformation efficiency	>1 x 10 ⁷ cfu/μg supercoiled DNA
Form	20 separate one-shot reactions
Biosafety level	S1 (Germany), BSL-2 (USA)
Stability	12 months after shipping
Storage	Store in cryo storage system at -90 °C to -60 °C
Shipping	Dry ice
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Application	Thaw a vial of competent <i>E. coli</i> TOP10 on ice. Add up to 10 μl DNA (e.g., from a StarGate [®] ligation reaction) to the thawed competent <i>E. coli</i> TOP10 cells. Mix gently (do not vortex) and then incubate for 30 min on ice. Repeat gentle mixing and incubate for 5 min at 37 °C. Mix gently and then incubate for 2-5 min on ice. Add 900 μl LB medium and shake for 45 min at 37 °C. This incubation step is necessary especially when using kanamycin to express resistance genes prior to plating on plates for selection. Plate 100 μl on LB agar containing antibiotic (if required) and 50 mg/L X-gal (optional). Centrifuge the residual 900 μl cell mixture for 30 sec in a microfuge, resuspend the cells with 100 μl LB medium and plate the whole amount as above. Incubate plates over night at 37 °C. Pick single colonies for further analyses (plasmid isolation, PCR...).

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

