

Data Sheet

Strep-Tactin®XT 4Flow® high capacity resin

50% suspension

Cat. No.: 2-5030-002, 2-5030-010, 2-5030-025,
 2-5030-500

Version: 5.0
 Revision Date: 09.01.2025

Description	Strep-Tactin®XT 4Flow® high capacity is a protein purification resin with a higher density of immobilized Strep-Tactin®XT compared to Strep-Tactin®XT 4Flow®. Strep-Tactin®XT is a streptavidin variant with optimized binding properties for Twin-Strep-tag® and Strep-tag®II fusion proteins.
Support	4% agarose, crosslinked
Form	50% suspension in buffer, pH 8.0 (100 mM Tris/HCl pH 8.0, 150 mM NaCl, 1 mM EDTA, 0.02% sodium azide)
Bead size	50-150 µm, spherical
Binding capacity	31 mg protein/ml resin Binding capacity was determined with a 50 kDa Twin-Strep-tag® fusion protein. Please note that the binding capacity is protein dependent.
Stability	12 months after shipping
Storage	2-8 °C
Shipping	room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.

Application	The suspension can be used for packing of gravity flow columns since the resin is optimized for column affinity chromatography, rather than for batch purification. To pack a column with 1 ml resin bed volume, 2 ml 50% suspension are necessary. Elution: 10x Buffer BXT (Cat. No. 2-1042-025, 1 M Tris/HCl pH 8.0, 1.5 M NaCl, 10 mM EDTA, 500 mM biotin) Regeneration: Buffer XT-R (Cat. No. 2-1045-250, 3 M MgCl ₂ , ready-to-use) or 100 mM NaOH
--------------------	--

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.