

Cell isolation with Strep-Tactin® Magnetic Microbeads

for PBMCs or other single cell suspensions

1. Required material

Cat. No.	Product	Required/total cells		
		1 x 10 ⁷	1 x 10 ⁸	1 x 10 ⁹
6-5510-050	Strep-Tactin® Magnetic Microbeads, 750 µl	15 µl	150 µl	1500 µl
-	Protein of choice fused to a Strep-tag®II or Twin-Strep-tag®	1 µg*	10 µg*	100 µg*
6-6325-001	Biotin stock solution, 100 mM, 1 ml	100 µl	150 µl	600 µl
-	Buffer for cell isolation, e.g. 1x PBS containing 1 mM EDTA and 0.5% BSA	30-40 ml	60-70 ml	180-190 ml

*Amounts refer to a 50 kDa Twin-Strep-tag® protein

2. Initial preparations

2.1. Reagent preparation

Volumes are suitable for isolating target cells out of up to 1 x 10⁷ PBMCs. For higher cell numbers, protein and Strep-Tactin® Magnetic Microbead volumes should be upscaled linearly according to total cell numbers (e.g., for 5 x 10⁷ cells use 5x indicated volume of Twin-Strep-tag® protein). Adapt other volumes according to Table 1.



Watch this How-to video to see an exemplary isolation:

<https://www.youtube.com/watch?v=Er-HVYXluH8>



Cell labeling and isolation (3.1. and 3.2.) has to be performed **at 4°C**. Please make sure that all reagents and cells are accordingly refrigerated before starting the protocol. The subsequent **removal of reagents and washing** (3.3 and 3.4) has to be performed **at room temperature**.

- 2.1.1.** Prepare a cell isolation buffer of choice, e.g. 1x PBS containing 1 mM EDTA and 0.5% BSA. The buffer should not contain biotin.
- 2.1.2.** *Optional:* Wash Strep-Tactin® Magnetic Microbeads before use to remove sodium azide. Add **1 ml** cell isolation buffer to required volume of microbeads (see 2.1.4.). Mix carefully and separate beads from buffer using a magnet. Discard supernatant and resuspend magnetic microbeads in cell isolation buffer (initial volume as in 2.1.4.).
- 2.1.3.** Dilute your Strep-tag®II or Twin-Strep-tag® fusion protein to a concentration of **50 µg/ml** with cell isolation buffer.



Titration of optimal isolation conditions might be necessary. The following instructions are an example for isolating cells with a 50 kDa protein fused to a Twin-Strep-tag®.

- 2.1.4.** Mix **20 µl** Twin-Strep-tag® fusion protein (2.1.3.) with **15 µl** Strep-Tactin® Magnetic Microbeads (vortex before pipetting!) in an Eppendorf tube. Incubate under constant gentle agitation for **5 min** (up to 24 h) **at 4 °C**.
- 2.1.5.** Prepare **1 mM** Biotin Elution Buffer by diluting **100 µl** of 100 mM Biotin stock solution in **10 ml** cell isolation buffer. Mix thoroughly. Keep at **room temperature**.

2.2. Sample preparation

Prepare 1×10^7 PBMCs in **30 µl** cell isolation buffer. Volume of cell isolation buffer should be upscaled linearly for higher cell numbers (e.g., use 5x 30 µl cell isolation buffer for 5×10^7 total cells). Cells should be cooled down to **4 °C** before starting the protocol.

Table 1: Recommended volumes & tube sizes for different cell numbers

Starting cell number	Recommended tube size [ml]	Resuspension volume [ml]	Total Biotin Elution Buffer [ml]	3.3.1. [ml]
$\leq 1 \times 10^7$	15	5	10	5
$\leq 1 \times 10^8$	15	10	15	7.5
$\leq 1 \times 10^9$	50	30	60	30

3. Protocol

3.1. Cell labeling



Perform all steps at **4 °C**.



3.1.1. Add PBMCs to the pre-incubated protein-Microbead preparation (2.1.4.) and mix thoroughly by gentle pipetting.



3.1.2. Incubate for **10 min** under gentle constant agitation, e.g. on a roller mixer, to prevent cells from sedimentation. Continue with 3.2.1.

3.2. Magnetic cell isolation




Perform all steps at **4 °C**.



3.2.1. Add the incubated cells to **5 ml** cell isolation buffer. Mix thoroughly by gentle pipetting.



3.2.2. Incubate the tube on a magnet for **1 min**, remove entire supernatant carefully.



3.2.3. Resuspend cells in **5 ml** cell isolation buffer. Incubate the tube on a magnet for **1 min** (see 3.2.2.), remove entire supernatant carefully.



3.2.4. Repeat step 3.2.3. once. Continue with step 3.3.1.

3.3. Removal of magnetic microbeads



Perform all steps at **room temperature**.



3.3.1 Resuspend cells in **5 ml** Biotin Elution Buffer (2.1.5.). Mix by thoroughly by pipetting and incubate for **5 min** at **room temperature** on a roller mixer.



3.3.2. Incubate the tube on a magnet for **1 min**, collect entire supernatant carefully and transfer it to a new collection tube.



3.3.3. Repeat steps 3.3.1. and 3.3.2. once.

3.3.4. Pool the supernatants and collect cells by centrifugation (400 x g, 6 - 10 min).

Remove supernatant and resuspend cells in appropriate buffer or medium for further applications.



If further removal of magnetic microbeads is needed (e.g. for further positive enrichment steps), proceed to step **3.4**.

3.4. Removal of remaining magnetic microbeads



Perform all steps at **room temperature**.



3.4.1 Discard supernatant carefully. Resuspend cell pellet in **5 ml** cell isolation buffer and incubate for **5 min** under agitation (e.g. on a roller mixer) at **room temperature**.



3.4.2. Place tube back on the magnet (to remove any potential residual beads) and incubate for **3 min**.

3.4.3. After incubation, transfer supernatant to a new tube and centrifuge cells for **6 – 10 min** at **400 x g**.

3.4.4. Remove supernatant and resuspend cells in appropriate buffer or medium for further applications.



To perform further positive isolation or depletion steps, please start the protocol once more at **3.1**.

4. Troubleshooting


Low yield	
<i>Option 1</i>	Titrate the ratio between Twin-Strep-tag® fusion protein and Strep-Tactin® Magnetic Microbeads for different cell numbers.
<i>Option 2</i>	Increase incubation time of cells with protein-Microbead mix (3.1.2.).
<i>Option 3</i>	Make sure that you carefully remove supernatants during incubation on the magnet (3.2.) without disrupting the binding of the microbeads to the magnet.
<i>Option 4</i>	Check for biotin contamination in your samples
Low purity	
Increase number of washing steps (3.2.)	
Microbeads contamination	
Make sure that you carefully remove supernatants during incubation on the magnet (3.3. and 3.4.) without disrupting the binding of the microbeads to the magnet.	
High amount of cell death	
Make sure that you always work at the recommended temperatures	

5. Example data

CD3⁺ T cell isolation

Isolation of CD3⁺ T cells from PBMCs using a Fab Fragment against CD3 fused to a Twin-Strep-tag®. Unlysed cells were stained with CD3-FITC (OKT-3) / CD45-APC (2D1) and analyzed by flow cytometry (CytoFlex, BC). Dead cells were excluded from the analysis using DAPI staining. Doublet and debris discrimination were performed using different FSC/SSC signals.

Before isolation



After isolation

