

Expression Vector (Acceptor Vector)

Host	<i>E. coli</i>		Mammalia			Yeast	Baculo	Name	Secretion	N-term	GOI	C-term	Resistance
Promoter	tet	T7	CMV	CMV	CMV	CUP1	Polyhedrin						
Vector	pASG	pPSG	pDSG (MEXi)	pCSG	pESG	pYSG	pLSG						
Cat.no.	5-4000-001	5-4200-001	5-5200-001	5-5000-001	5-4400-001	5-4600-001	5-4800-001	IBAwt1	no				Amp
	5-4005-001	5-4205-001		5-5005-001	5-4405-001	5-4605-001	5-4805-001	IBA5	no				Amp
	5-4105-001	5-4305-001	5-5222-001	5-5105-001	5-4505-001	5-4705-001	5-4905-001	IBA105	no				Amp
	5-4035-001	5-4235-001		5-5035-001	5-4435-001	5-4635-001	5-4835-001	IBA35	no				Amp
	5-4025-001	5-4225-001		5-5025-001	5-4425-001	5-4625-001	5-4825-001	IBA25	no				Amp
	5-4003-001	5-4203-001		5-5003-001	5-4403-001	5-4603-001	5-4803-001	IBA3	no				Amp
	5-4103-001	5-4303-001	5-5220-001	5-5103-001	5-4503-001	5-4703-001	5-4903-001	IBA103	no				Amp
	5-4033-001	5-4233-001	5-5209-001	5-5033-001	5-4433-001	5-4633-001	5-4833-001	IBA33	no				Amp
	5-4045-001	5-4245-001	5-5214-001	5-5045-001	5-4445-001	5-4645-001	5-4845-001	IBA45	no				Amp
	5-4145-001	5-4345-001		5-5145-001	5-4545-001	5-4745-001	5-4945-001	IBA145	no				Amp
	5-4043-001	5-4243-001	5-5211-001	5-5043-001	5-4443-001	5-4643-001	5-4843-001	IBA43	no				Amp
	5-4143-001	5-4343-001		5-5143-001	5-4543-001	5-4743-001	5-4943-001	IBA143	no				Amp
	5-4023-001	5-4223-001		5-5023-001	5-4423-001	5-4623-001	5-4823-001	IBA23	no				Amp
	5-4123-001	5-4323-001		5-5123-001	5-4523-001	5-4723-001	5-4923-001	IBA123	no				Amp
	5-4065-001	5-4265-001		5-5065-001	5-4465-001	5-4665-001	5-4865-001	IBA65	no				Amp
	5-4063-001	5-4263-001		5-5063-001	5-4463-001	5-4663-001	5-4863-001	IBA63	no				Amp
	5-4062-001	5-4262-001		5-5062-001	5-4462-001	5-4662-001	5-4862-001	IBA62	no				Amp
	5-4162-001	5-4362-001		5-5162-001	5-4562-001	5-4762-001	5-4962-001	IBA162	no				Amp
	5-4064-001	5-4264-001		5-5064-001	5-4464-001	5-4664-001	5-4864-001	IBA64	no				Amp
	5-4164-001	5-4364-001		5-5164-001	5-4564-001	5-4764-001	5-4964-001	IBA164	no				Amp
	5-4167-001	5-4367-001		5-5167-001	5-4567-001	5-4767-001	5-4967-001	IBA167	no				Amp
	5-4168-001	5-4368-001		5-5168-001	5-4568-001	5-4768-001	5-4968-001	IBA168	no				Amp
	5-4001-001		5-5201-001	5-5001-001	5-4401-001		5-4801-001	IBAwt2	yes				Amp
	5-4004-001							IBA4	yes				Amp
	5-4104-001		5-5221-001	5-5104-001	5-4504-001		5-4904-001	IBA104	yes				Amp
	5-4002-001							IBA2	yes				Amp
	5-4102-001		5-5219-001	5-5102-001	5-4502-001		5-4902-001	IBA102	yes				Amp
	5-4044-001							IBA44	yes				Amp
	5-4144-001			5-5144-001	5-4544-001		5-4944-001	IBA144	yes				Amp
				5-5142-001	5-4542-001		5-4942-001	IBA142	yes				Amp

Secretion signal in *E. coli*: OmpA
 Secretion signal in eukaryotes: BM40

