

## Data Sheet

# 10x Buffer W Strep-Tactin®/XT Wash Buffer

Cat. No.: 2-1003-100

Version: 5.2

Lot No.: 1003-

Revision Date: 25.02.2020

<b>Description</b>	Buffer W is a Strep-Tactin®/XT Wash Buffer for the preparation of protein extracts and wash of immobilized Strep-Tactin® and Strep-Tactin®XT. Buffer is supplied in 10x concentration.
<b>Components (10 x concentrated solution)</b>	<b>10x Buffer W</b> (Strep-Tactin®/XT Wash Buffer), pH 8,0 1 M Tris-Cl 1.5 M NaCl 10 mM EDTA
<b>Stability</b>	12 months after shipping
<b>Storage</b>	2-8 °C for both concentrated and working solution
<b>Shipping</b>	room temperature
<b>Hazards</b>	The concentrated buffer is classified as hazardous according to (EC) No 1272/2008 [CLP]. See provided Material Safety Data Sheet for details. The working solution is not hazardous.
<b>Working solution</b>	Dilute 1 volume of concentrated solution (10x) with 9 volumes of demineralized or distilled water.



Go digital and help the environment. Please download all up-to-date manuals, protocols and other material from <http://www.iba-lifesciences.com>.

### For research use only

#### Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.