



IBA Lifesciences GmbH
Rudolf-Wissell-Str. 28
37079 Goettingen
Germany
Tel.: +49 (0) 551-5 06 72-0
E-mail: info@iba-lifesciences.com
www.iba-lifesciences.com

Data Sheet

GFP-Twin-Strep-tag control protein

Cat. No.: 2-1007-005

Version: 3.1
Revision Date: 28.07.2021

Description	GFP-Twin-Strep-tag is a control protein consisting of the green fluorescent protein (GFP) C-terminally tagged with the Twin-Strep-tag.
MW	30 kDa
Form	Dissolved in 1x Buffer W
Concentration	1 mg/ml
Amount	500 µl
Stability	12 months after shipping
Storage	-25 °C to -15 °C
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Application	It is suitable for test purifications, detections, and immobilizations with Strep-Tactin® or Strep-Tactin®XT.

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.