



**IBA Lifesciences GmbH**  
Rudolf-Wissell-Str. 28  
37079 Goettingen  
Germany  
Tel.: +49 (0) 551-5 06 72-0  
E-mail: [info@iba-lifesciences.com](mailto:info@iba-lifesciences.com)  
[www.iba-lifesciences.com](http://www.iba-lifesciences.com)

## Data Sheet

# Strep-Tactin® PE

Cat. No.: 6-5000-001

Version: 7.4  
Revision Date: 29.09.2021

<b>Description</b>	Strep-Tactin® conjugated to R-Phycoerythrin (PE) suitable for cell staining.
<b>Form</b>	Dissolved in DPBS (Dulbecco's phosphate buffered saline) without Ca & Mg, contains KCl; pH 7.0 – 7.5
<b>Amount</b>	50 µl
<b>Stability</b>	12 months after shipping; handle under <b>sterile</b> conditions.
<b>Storage</b>	2-8 °C; <b>Do not freeze!</b>
<b>Shipping</b>	Blue ice
<b>Hazards</b>	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
<b>Application</b>	Strep-Tactin® fluorescently conjugated to either R-Phycoerythrin (PE) or Allophycocyanin (APC) is suitable for staining or fluorescent sorting of antigen-specific cytotoxic T cells (MHC I Streptamer®) or surface marker-specific cells (Fab-TACS®/Nano-TACS®) based on our Strep-tag® technology. For this purpose, fluorescent Strep-Tactin® has to be combined with Twin-Strep-tag®-containing MHC I complexes (MHC I-Streps) or low affinity Fab fragments/nanobodies (Fab-Streps/Nano-Streps). Biotin causes the dissociation of all reagents from the cells, yielding label-free populations for unbiased further use.
<b>Recommended usage</b>	1 µl per 2 x 10 <sup>6</sup> (MHC I-Streps) or 5 x 10 <sup>6</sup> (Fab-Streps) total cells

### For research use only

#### Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

#### Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at [info@iba-lifesciences.com](mailto:info@iba-lifesciences.com) or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.