



IBA Lifesciences GmbH
Rudolf-Wissell-Str. 28
37079 Goettingen
Germany
Tel.: +49 (0) 551-5 06 72-0
E-mail: info@iba-lifesciences.com
www.iba-lifesciences.com

Data Sheet

Strep-Tactin® AP

Cat. No.: 2-1503-001

Version: 8.0
Revision Date: 15.06.2021

Description	Strep-Tactin® labelled with alkaline phosphatase (AP, min. 2 moles per mole Strep-Tactin®). Strep-Tactin®* is a streptavidin variant specially optimized for Strep-tag®II and Twin-Strep-tag® fusion protein binding.
Form	Suspension in buffer (50 mM Tris-HCl, 0.1 M NaCl, 1 mM MgSO ₄ , 0.1 mM ZnCl ₂ , 0.1% (w/v) NaN ₃ , 10 mg/ml bovine serum albumin, pH 8)
Alkaline phosphatase activity	1500-2000 units/mg (DEA units, unit definition: one unit is the enzyme activity which dephosphorylates 1 µmole of orthophosphoric monoester per min under assay conditions (25 °C, pH 9.8))
Amount	500 µl
Concentration	0.175-0.185 mg/ml conjugate (Bradford dye binding assay)
Stability	12 months after shipping
Storage	Recommended: 2-8 °C
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Application	Detection on microplates: dilute 1:800 Western blot: dilute 1:4000 or higher for chromogenic detection

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.