



**IBA Lifesciences GmbH**  
Rudolf-Wissell-Str. 28  
37079 Goettingen  
Germany  
Tel.: +49 (0) 551-5 06 72-0  
E-mail: [info@iba-lifesciences.com](mailto:info@iba-lifesciences.com)  
[www.iba-lifesciences.com](http://www.iba-lifesciences.com)

## Data Sheet

# Strep-tag®II peptide

Cat. No.: 2-1018-002

Version: 5.1  
Revision Date: 28.07.2021

<b>Description</b>	<b>Synthetic Strep-tag®II peptide recognized by Strep-Tactin®, Strep-Tactin®XT as well as the monoclonal antibodies, StrepMAB-Classic and StrepMAB-Immo.</b>
<b>Sequence</b>	SAWSHPQFEK
<b>MW</b>	1216 Da
<b>Form</b>	Lyophilized as TFA salt.
<b>Amount</b>	1.8 mg
<b>Stability</b>	6 months after shipping
<b>Storage</b>	Recommended: -25 to -15 °C
<b>Shipping</b>	Room temperature
<b>Hazards</b>	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
<b>Application</b>	Dissolve the whole amount in 3 ml of a suitable buffer to get a 0.5 mM solution. It is recommended to use the same buffer as for washing the column. It can be used for competitive elution of Strep-tag®II fusion proteins or as control protein in dot blot assays.

### For research use only

#### Trademark information

The owners of trademarks marked by “®” or “™” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

#### Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at [info@iba-lifesciences.com](mailto:info@iba-lifesciences.com) or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.