



IBA Lifesciences GmbH
Rudolf-Wissell-Str. 28
37079 Goettingen
Germany
Tel.: +49 (0) 551-5 06 72-0
E-mail: info@iba-lifesciences.com
www.iba-lifesciences.com

Data Sheet

10x Buffer CI

Cat. No.: 6-6320-085

Version: 1.2
Revision Date: 04.10.2021

Description	10x PBS containing 10 mM EDTA and 5% BSA, sterile filtered.
Amount	85 ml
Stability	12 months after shipping; handle under sterile conditions.
Storage	2-8 °C; may also be stored at 15-25 °C
Shipping	Room temperature or blue ice
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Application	10x Buffer CI is suitable for our Strep-tag® technology approaches for cell isolation and staining (Fab-TACS®/Nano-TACS® and MHC I Streptamer®).
Usage	Dilute 1:10 with ddH ₂ O for 1x Buffer CI working solution.

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.