



IBA Lifesciences GmbH
Rudolf-Wissell-Str. 28
37079 Goettingen
Germany
Tel.: +49 (0) 551-5 06 72-0
E-mail: info@iba-lifesciences.com
www.iba-lifesciences.com

Data Sheet

10x Buffer CI

Cat. No.: 6-6320-025, 6-6320-085

Version: 1.3
Revision Date: 28.02.2024

| | |
|--------------------|---|
| Description | 10x PBS containing 10 mM EDTA and 5% BSA, sterile filtered. |
| Stability | 12 months after shipping; handle under sterile conditions. |
| Storage | 2-8 °C; may also be stored at 15-25 °C |
| Shipping | Room temperature or blue ice |
| Hazards | Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided. |
| Application | 10x Buffer CI is suitable for our Strep-tag® technology approaches for cell isolation and staining. |
| Usage | Dilute 1:10 with ddH ₂ O for 1x Buffer CI working solution. |

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.